

Seminal plasma as a noninvasive sensitive diagnostic matrix for the evaluation of PFAS exposure in firefighters

Jana Navratilova¹, Michal Jeřeta², Ales Pindur³, Petr Šenk¹, Pavel Čupr¹

¹RECETOX, Faculty of Science, Masaryk University, Kamenice 753/5, 625 00 Brno, Czech Republic, Faculty of Sports Studies, Masaryk University, Kamenice 753/5, 625 00 Brno, Czech Republic

²Center of Assisted Reproduction, Department of Gynecology and Obstetrics, Faculty of Medicine, Masaryk University Brno and University Hospital Brno, Czech Republic

³Training Center of Fire Rescue Service, General Directorate of Fire Rescue Service of the Czech Republic, Ministry of the Interior, Trnkova 85, 628 00 Brno, Czech Republic

Introduction


- Per- and polyfluoroalkyl substances (PFASs) are a group of synthetic chemicals with water resistant properties primarily used in variety of products such as carpets, furniture, paints, fabrics etc.
- PFAS exposures have been associated with adverse health outcomes, including cancer and fertility issues
- Firefighters have elevated concentrations of certain PFASs compare to non-firefighters
- Firefighters are exposed to PFAS through the combustion of PFAS-containing products such as furniture and carpet and through firefighting gear and firefighting foams that contain these compounds
- Seminal plasma is an alternative matrix to blood to monitor for PFAS exposure

Aims and Objectives

- To optimize sample preparation method for selective extraction of 15 PFASs from seminal fluid
- To develop and validate a method for simultaneous quantification of PFASs using liquid-chromatography tandem mass spectrometry (LC-MS/MS)

Materials and methods

- A 96 well plate sample preparation method:**



- *1% formic acid solution in acetonitrile
- **The eluent was evaporated to dryness under nitrogen at 50°C and reconstituted with 100µL of water-methanol (50:50 v/v)

Figure 1: Workflow for detection of PFASs based on a 96-well plate sample preparation and analysis on LC-ESI-MS/MS.

- LC-MSMS method:**

➤ Chromatographic conditions

Separation column: SYNERGI 4µ Fusion Max-RP 80Å 100 mm x 2 mm (Phenomenex, USA)

Elution conditions: Gradient elution with a mixture of 5mM ammonium acetate and methanol; 400 µl/min flow rate

➤ MS/MS detection

QTrap 5500 (ABSciex, CA, USA) operating in ESI-negative mode

Results

- The method validation was performed based on FDA guidelines
- Validation criteria included limit of detection (LOD), limit of quantitation (LOQ), accuracy and precision
- Due to the lack of the certified reference standard material the validation criteria of the method were assessed in pool of in house spiked seminal samples (QC)




Figure 2: Precision and accuracy* of the extraction and quantification method

* The accuracy of the sample preparation method was evaluated as a percent recovery of the amount of target analyte added into the sample




Figure 3: Representative LCMSMS chromatogram of PFAS mix calibration standard 1ng/mL

Figure 4: LCMSMS chromatogram of real sample prepared on 96 well plate. The x-axis is time (min) from 0.0 to 4.0. The y-axis is intensity from 0.0 to 9.0e4. Multiple peaks are visible, with the most prominent at 0.67 min.

Conclusions

- ✓ Selective sample preparation method was optimized for detection of PFAS in seminal fluid samples by LC-MS/MS
- ✓ The precision and accuracy of the in house prepared QC samples showed that the method is robust and accurate
- ✓ Capability of being able to quantitatively measure PFASs in seminal fluid in pg/mL concentration range can help to identify and establish their role in male infertility and serve as alternative method for their monitoring

Acknowledgement

This project has received funding from the European Union's Horizon 2020 research and innovation programme under grant agreement No 101003355 and grant agreement No 857340. The project was supported by MH CZ – DRO (FNBr; 65269705) and project AZV NV18-01-00544. Authors also thank to Research Infrastructure RECETOX RI (No LM2018121) financed by the Ministry of Education, Youth and Sports, and Operational Programme Research, Development and Innovation - project CETOCOEN EXCELLENCE (No CZ.02.1.01/0.0/0.0/17_043/0009632) for supportive background. This presentation reflects only the author's view, and the European Commission is not responsible for any use that may be made of the information it contains.